

Push-Pull: KDS 120



Tomopal Part #	Manufacturer #	Voltage
781-1203 USA	781120 USA	115V~, 0.1A
789-1203 CE	789120 CE	230V~, 0.05A

Features:

- Two Syringe Push-Pull Pump One + One 10 µl-15 ml each
- This pump provides simultaneous infusion and withdrawal at the same rate with opposing syringes on the same drive screw. The push-pull mode is designed for one cycle only.
- Dispense Volume Mode: The pump keeps track of volume dispensed and automatically stops when target volume is reached.
- Run Mode: Pump runs at the set flow rate until stopped manually.
- Min. Flow Rate (10 µl syringe): 0.1 µl/hr
- Max. Flow Rate (10 ml syringe): 127 ml/hr