

Infusion/Withdrawal Pumps: KDS 210



Tomopal Part #	Manufactuerv #	Voltage
781-2103 USA	781210 USA	115V~, 0.1A
789-2103 CE	789210 CE	230V~,0.05A

Features:

- Two-Syringe Infusion/Withdrawal Pump 10 µl-140 ml each
- A target volume can be entered for infusion and refill independently, and the pump automatically stops when this volume is reached. Pump displays an initial volume of zero and increases as the dispense proceeds to the target volume. The target volume can be reviewed or changed as the pump continues to operate.
- RS-232 Interface
- TTL Interface
- Min. Flow Rate (10 µl syringe): 0.001 µl/hr
- Max. Flow Rate
 - 10 ml Syringe: 21 ml/min
 - 60 ml Syringe: 86 ml/min
 - 140 ml Syringe: 145 ml/min